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These instructions detail the care and use of the GreenLight HPS fiber optic. They are not recommendations for the medical or surgical application of this device. Any physician using this device should be thoroughly familiar with the surgical procedure being performed prior to using the device. This Fiber is designed for use with compatible AMS GreenLight laser systems. Refer to the laser system Operator Manual, Safety and Professional Information Sections, for specific instructions concerning warnings, cautions, contraindications and clinical use of the laser.

Package Contents:

Description **REF** 10-2090

(1) Individually packaged, sterile, disposable GreenLight HPS fiber optic with fiber card

Description:

The GreenLight HPS fiber optic is a sterile, single-use, disposable delivery system with a 600 micron fiber. The device is designed to deliver laser energy at the 532nm wavelength. The proprietary connector permits the laser to recognize the delivery device and activates the software safety features.

The GreenLight HPS fiber optic is designed to deliver laser energy from the tip laterally at 70 degrees. The GreenLight HPS fiber optic is intended to be used in surgical applications where lateral delivery of the laser energy is needed (i.e. procedures where lateral delivery of laser energy is normally accomplished with fibers alone or in conjunction with mirrors or fiber diverters).

The GreenLight HPS fiber optic is designed to accommodate a variety of surgical applications and endoscopic instrumentation. The GreenLight HPS fiber optic can be used in contact or non-contact and rotated 360 degrees which allows tissue access in multiple planes and is compatible with standard endoscopes, cystoscopes and Albarran bridges. This device will cut, vaporize, or coagulate tissue in liquid endoscopic surgical procedures.

Instructions for Use:

CAUTION
Check fiber package for damage. If damage is noted or sterility is in question, DO NOT USE.

- Open box remove fiber container. Remove GreenLight HPS fiber optic from its sterile container using aseptic technique. Bend back package tabs fully prior to removing fiber from packaging nest.
- Turn laser system on. Wait for self-test to be completed and the prompt "Insert Fiber Card" to appear.
- Remove fiber card from pouch and insert into the laser console (see Operators Manual). The image on the card should be facing the front of the laser. After the card is read, there will be a prompt to "Attach Device."

CAUTION
The fiber card must remain inserted until the procedure is complete. Premature removal of the card could result in the discontinued use of the fiber during the procedure.

- Insert the GreenLight HPS fiber optic into the fiber port on the laser system and turn it ¼ turn clockwise, until it locks.
- Before continuing with the surgical procedure, the GreenLight HPS fiber optic should be checked for damage during shipment. Place the laser in READY mode to activate the aim beam.

CAUTION
Do not press the footswitch while checking the aim beam.

Place the distal end of the device on a non-reflective sterile surface and turn slightly until the aim beam can be visualized. If the aim beam is not seen, the device may be defective and should not be used. Once the GreenLight HPS fiber optic has been checked, return the laser to the STANDBY mode.

- Position the GreenLight HPS fiber optic at the selected treatment site. If an endoscope is used, the tip of the GreenLight HPS fiber optic should be in clear view and extended beyond the end of the endoscope. The output beam of the GreenLight HPS fiber optic is aligned with the raised surface on the circular knob, therefore, the direction of the laser energy can always be determined.
- Place the laser in READY mode to enable the footswitch control.
- Do not exceed the following maximum power recommendations:

Environment	Wavelength	Max. Power
Liquid	532 nm	120 watts

The recommended maximum powers shown above are based on results obtained during clinical evaluations. They are not necessarily indicative of clinical performance. Exceeding these powers may result in product damage.

- AMS software automatically corrects for fiber losses such that the power level shown on the screen video display indicates the amount of power delivered to tissue.
- Treatment times vary based on distance to tissue, power settings, and other factors.
- If, during the procedure, the tip accumulates debris, turn the laser to the STANDBY mode and carefully wipe the tip clean with a wet sterile gauze or towel. Wipe the fiber starting at the tip and moving down the fiber.

WARNING
Never activate the laser energy if the fiber tip does not extend beyond the end of the endoscope.

- The fiber can be used up to an energy or time limit. The user is alerted by an audible tone and a screen message when less then 50,000 joules or 5 minutes remain. To continue a procedure when a limit has been reached, a new fiber card and fiber must be inserted. *Note: Do not cleave or trip the GreenLight HPS fiber optic.*
- If excessive tissue or debris is allowed to accumulate on the fiber cap, over heating of fiber cap will occur leading to premature fiber degradation and/or fiber failure.
- Dispose of the GreenLight HPS fiber optic following Vacho or OSHA protocol.

PROCEDURAL WARNINGS
To avoid GreenLight HPS fiber optic damage or breakage:
- Do not bury tip in tissue. - Do not retract or probe tissue with the device. - Do not bend at sharp angles. - If using a deflecting mechanism, such as an Albarran bridge or a fiber diverter, do not attempt to remove the GreenLight HPS fiber optic while the bridge is deflected. - In the unlikely event of a detached tip, it may be visually located through an appropriate scope and removed using forceps. Irrigate the area thoroughly to remove any traces of fiber or other material. - Reuse or misuse of a GreenLight HPS fiber optic will result in an increased potential for damage to the GreenLight HPS fiber optic and ancillary equipment, i.e., cystoscopes.
USAGE WARNINGS
- A sudden increase in bubbles at the probe tip indicates insufficient irrigation, tissue contamination or device failure. If increased irrigation does not eliminate the bubbling then stop, clean and inspect the device. Discontinue use and discard instrument if damage is evident or bubbling occurs upon - Reuse of a GreenLight HPS fiber optic can pose health and safety risks to patients and may compromise treatment outcomes. - Do not use a damaged fiber optic.
FIRE HAZARD WARNINGS
- The GreenLight HPS fiber optic is a precision device. Damage to the fiber can result in the emission of uncontrolled laser energy. To prevent draps fires or burns: - Do not wrap any portion of the fiber in drapes or cloth. - Do not attach the fiber directly to surgical draps with a crushing clamp such as a hemostat. - Do not place or drop instrumentation onto the fiber.

How Supplied And Storage

WARNING: Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, call your AMS representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or reesterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or reesterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Made in USA

GreenLight HPS® is a trademark of American Medical Systems.

P/N 0127-1410 Rev. C

Warranty:

The Warranty for this product covers only defects in the materials and workmanship. The GreenLight HPS fiber optic is a disposable device and is intended to be used in a single procedure and then discarded. Should you need to return this product, call the AMS Customer Care Center or the appropriate international distributor for "Return Authorization." Instructions for returning product to AMS are detailed in the Accessory Return/Repair Policy form included with each device.

Accessory Return/Repair Policy
Return/Repair Policy and Procedures for End-Users of AMS Accessories

Overview of AMS's Return/Repair Policy:

Generally, any AMS accessory that fails during its warranty period because of defects in materials or workmanship—NOT MISUSE—may be returned by the customer to AMS for replacement. Additionally, certain accessories that fail after the warranty period has expired may be returned for repair. Warranty policy for returned products may be obtained from your Laser System Operator's Manual. After warranty, the customer may be responsible for repair costs.

The policies and procedures outlined below are not a complete explanation of AMS's return and repair policies and procedures. Please refer to the AMS Terms and Conditions of Sale (on the back of the AMS Sales Order) for a complete explanation or contact the AMS Customer Care Center for this information.

Procedure For Returning AMS Accessories:

All products returning to AMS must have a signed Sterilization Certificate and a Return Materials Authorization (RMA) number. The Sterilization Certificate must be attached to the outside of the returning package. RMA numbers and Sterilization Certificates may only be obtained by calling the Customer Care Center at 1-800-356-7000 Monday through Friday, 7:00 am to 4:30 pm Pacific Standard Time. If you return a failed product that is under warranty, AMS will credit your account or replace the product only if you first obtain an RMA number. Products not covered by warranty must also have an RMA number, and, upon receipt, your return request will go through AMS's determination process before being accepted for credit consideration.

If a Warranty Return Form or any return questionnaire was supplied with the device or is provided by AMS's Customer Care Center when advised of the intent to return, the form must be completed and returned with the product.

Sterile/Non-Sterile Products:

In order to protect everyone who may come in contact with a "Returned Product," a Sterilization Certificate must be attached to the outside of the returning package.

If a product has not been used and is in its original unopened package, it may be returned non-sterile once an RMA number has been issued. Products must be returned within 30 days of the original shipment date to be given credit consideration. A 20% restocking fee will be charged for all unused and unopened returned products.

Products returned without the Sterilization Certificate on the outside of the returning package will be returned unopened. Boxes that are opened and found to be inappropriately decontaminated shall be contained and isolated from all other materials. The customer shall be notified by Customer Care of charges levied for the disposal of goods in accordance with OSHA Standards on goods will be repackaged in accordance with DOT regulations and returned to the customer as hazardous material at customer's expense.

Unused and Unopened Disposable Accessories:

Sterile disposable products can only be returned if the seal on the original package is not broken. A 20% restocking fee will be charged for all unused and unopened disposable accessories. Products must be returned within 30 days of the original shipment date to be given credit consideration.

Obtaining a Return Material Authorization (RMA) Number:

When calling AMS's Customer Care Center to request an RMA number, please be prepared to provide the following information:

- Your facility: Customer name, address and contact person.
- Your AMS Customer Number.
- Your AMS Sales Order Number.
- Your original Purchase Order Number.
- Catalog Number and quantity of product to be returned.
- Serial or Lot Number of product to be returned.
- Reason for return (including laser system error codes displayed, if available).
- Your "replacement" Purchase Order Number.

Shipping Procedures:

Ship the package, postage pre-paid. Indicate the RMA number on the outside of the package. Insure the package for its appropriate value (as indicated on the original purchase order), as AMS will not be financially responsible for the loss of any product being returned.

Note: If the product is not received by AMS within 30 days of original shipment date, the return product will not be accepted for credit consideration.

Replacement Products Out of Stock or Not Timely Delivered:

If, as a result of unforeseen events, AMS does not have the appropriate replacement product in stock or the replacement product is not timely delivered, AMS assumes no responsibility for monetary loss or damage resulting to the customer/end-user.

Determining if Credit is Due:

Within 30 days of receipt of a returned product, AMS will evaluate and test the returned item to determine if credit is due to the customer. AMS warrants products only for defects in materials and workmanship and only for a specific period, so it is possible that returned products may not be credited to a customer's account.